

Dexdomitor Doses

DOGS

Premed for routine surgeries

Combine _ label dose with opioid of choice at standard dose and administer IM 15-20 minutes prior to induction. Induction drug volume should be reduced to _ usual amount or less. Inhalant anesthesia should also be reduced by 50% on average.

Protocol for short non or mildly painful procedures (Penn Hip)

Combine _ label dose Dexdomitor with 0.2mg/kg Torb and administer IM or IV

Can reverse with equal volume (to Dexdomitor) of Antisedan IM

Mini Micro Rescue dose for rough recovery

- Administer 0.05-1ug/kg IV. This will be about 0.1-0.2 ccs for a 50 lb dog. Provides about 30 minutes of sedation to transition smoothly from anesthesia.
- Can repeat as needed or deliver as a CRI at 1-3 ug/kg/hr

CATS

Kitty Magic for Surgical or painful procedures

In a 10 lb cat combine and administer IM:

0.2cc Dexdomitor

0.2cc Ketamine

0.2 cc Buprenorphine (or other opioid)

Provides 30 minutes of profound sedation and analgesia typically sufficient to perform castration or less painful procedures or intubation. Occasionally small amounts of inhalant anesthesia by mask are required.

For simple non painful procedures dose can be reduced to

0.1cc Dexdomitor per 10 lbs

0.1cc Ketamine per 10 lbs

0.1cc Torbugesic per 10 lbs

RABBITS

Bunny Magic combine in the same syringe and administer IM:

0.05 mg/kg Dexdomitor

0.05 mg/kg buprenorphine

5 mg/kg ketamine.

As with kitty magic, double the ketamine for deeper anesthesia.